SIR-Spheres® Microspheres (Yttrium-90 Microspheres) SIR-Y001

ENGLISH

1. DESCRIPTION

SIR-Spheres microspheres consist of biocompatible microspheres containing yttrium-90. Yttrium-90 is a high-energy pure beta-emitting isotope with no primary gamma emission. The maximum energy of the beta particles is 2.27 MeV with a mean of 0.93 MeV. The half-life is 64.1 hours. The maximum range of emissions in tissue is 11 mm with a mean of 2.5 mm.

SIR-Spheres microspheres are implanted into hepatic tumors by injection into the hepatic artery using a catheter. The SIR-Spheres microspheres distribute non-uniformly in the liver, primarily due to the unique physiological characteristics of the hepatic arterial flow, the tumor to normal liver ratio of the tissue vascularity, and the size of the tumor. The tumor usually gets higher density per unit distribution of SIR-Spheres microspheres than the normal liver. The density of SIR-Spheres microspheres in the tumor can be as high as 5 to 6 times of the normal liver tissue. In general, 1 GBq of yttrium 90/kg of tissue provides 49.67 Gy of radiation dose. Once SIR-Spheres microspheres are implanted into the liver, they are not metabolized or excreted, and they stay permanently in the liver. Each device is for single patient use.

2. INTENDED USE

SIR-Spheres Y-90 resin microspheres are intended for implantation into hepatic tumors via the hepatic artery.

3. INDICATIONS FOR USE

SIR-Spheres Y-90 resin microspheres are indicated for the treatment of unresectable hepatocellular carcinoma (HCC) or unresectable metastatic liver tumors from primary colorectal cancer in patients refractory to or intolerant of chemotherapy.

4. HOW SUPPLIED

SIR-Spheres microspheres are provided in a vial with water for injection. Each vial contains 3 GBq of $Y^{90} \pm 10\%$ (at the time of calibration) in a total of 5 cc water for injection. Each vial contains 40 – 80 million microspheres with a diameter between 20 and 60 micrometers. The vial is shipped within a 6.4 mm minimum thickness lead pot. The package consists of a crimp sealed SIR-Spheres microspheres glass vial within a lead pot, and a package insert within Type A package.

Patient specific activities are prepared in accordance with Section 11.

The vial and its contents should be stored inside its transportation container at room temperature ($15-25^{\circ}\text{C}$; $59-77^{\circ}\text{F}$).

The calibration date (for radioactive contents) and the expiration information are printed on the vial label. The useful life of the SIR-Spheres microspheres ends 24 hours after the time of calibration.

5. ACCESSORIES

Accessories that may be used for the dose preparation and implant procedure include:

- Delivery Set (SIR-D001)
- V-Vial and capped 21 G needle (SIR-V001)
- V-Vial Holder (SIR-H001)
- Delivery Box (SIR-B001)
- Syringe Shield (SIR-S001)

These accessories are designed to shield the user during dose preparation as well as the user and patient during the implant procedure. Sections 12 - 15 provide additional safety precautions for users to ensure exposure to radiation is kept as low as reasonably achievable (ALARA).

6. CONTRAINDICATIONS

SIR-Spheres microspheres are contraindicated in patients who have:

- had previous external beam radiation therapy to the liver
- ascites or are in clinical liver failure
- markedly abnormal liver function tests (LFTs)
- Total bilirubin > 2.0 mg/dL and/or Albumin < 3.0 g/dL
- ≥ 30 Gy absorbed radiation dose to the lungs, as estimated by the Technetium-99m macroaggregated-albumin (^{99m}Tc MAA) study described in Section 10
- pre-assessment angiogram that demonstrates potential pathways for deposition of microspheres to non-target organs such as the stomach, pancreas, or bowel
- Are pregnant

7. WARNINGS

• Non-Target Delivery of SIR-Spheres microspheres

Inadvertent delivery of SIR-Spheres microspheres to extra-hepatic structures such as the stomach, duodenum, gallbladder, or pancreas may result in radiation injury to these structures, including but not limited to acute abdominal pain, acute gastritis, acute cholecystitis, acute pancreatitis, and peptic ulceration. Angiographic technique must be employed to prevent the non-target delivery of SIR-Spheres microspheres to any extra-hepatic structures.

• Radioembolization Induced Liver Disease (REILD)

Delivery of excessive radiation to the normal liver parenchyma may result in REILD. The risk of REILD may also be increased in patients with pre-existing liver disease. Consideration should be given to reducing the prescribed activity of SIR-Spheres microspheres in the following clinical settings¹:

- Reduced liver functional reserve due to steatosis, steatohepatitis, hepatitis or cirrhosis
- Elevated baseline bilirubin level
- Prior hepatic resection
- Prior liver directed therapy
- Extensive prior treatment with systemic chemotherapy and/or biologic therapies

• Radiation Pneumonitis

High levels of implanted radiation and/or excessive shunting to the lung may lead to radiation pneumonitis. The lung radiation dose must be limited to ≤ 30 Gy for a single treatment session and ≤ 50 Gy cumulative dose.

• Other Tumor Types

The safety, performance, and benefit/risk profile of SIR-Spheres in treating particular types of tumors outside the indications for use has yet to be established.

8. PRECAUTIONS

- Safety and effectiveness of this device in pregnant women, nursing mothers or children have not been established.
- A SPECT or PET scan of the upper abdomen be performed immediately after implantation of SIR-Spheres microspheres. The SPECT or PET scan will detect radiation from the yttrium-90 to confirm placement of the microspheres in the liver.
- This product is radioactive. The use of this device is regulated, and local regulations must be followed when handling this device.
- Accepted radiation protection techniques should be used to protect staff when handling both the isotope and the patient.
- Patients may experience gastric problems following treatment, but proton pump inhibitors (PPI) or histamine H2-receptor antagonists (H-2 blocking agents) may be used the day before implantation of SIR-Spheres microspheres and continued as needed to reduce gastric complications.
- Patients may experience abdominal pain immediately after administration of SIR-Spheres microspheres and pain relief may be required.
- SIR-Spheres microspheres demonstrated a mild sensitization potential when tested dermally in an animal model.

9. ADVERSE EVENTS

When the patient is treated with proper technique, without excessive radiation to any organ, the common adverse events after receiving the SIR-Spheres microspheres are fever, transient decrease of hemoglobin, transient thrombocytopenia, mild to moderate abnormality of liver function tests (mild increase in aspartate aminotransferase, alkaline phosphatase, bilirubin), abdominal pain, nausea, vomiting, and diarrhea.

Potential Serious Adverse Events Due to High Radiation

Acute Pancreatitis: causes immediate severe abdominal pain. Verify by SPECT or PET of the abdomen and test for serum amylase.

Radiation Pneumonitis: causes excessive non-productive cough. Verify by X-ray evidence of pneumonitis.

Acute Gastritis: causes abdominal pain. Verify by standard methods to diagnose gastric ulceration.

Acute Cholecystitis: causes significant upper abdominal pain and may require cholecystectomy for resolution. Verify by appropriate imaging studies.

Radioembolization induced liver disease (REILD): REILD is a rare complication following Selective Internal Radiation Therapy (SIRT). REILD is characterized by a well-defined constellation of temporal, clinical, biochemical, and histopathologic findings. It typically manifests approximately 4 to 8 weeks post-SIRT and is characterized clinically by jaundice and ascites in the absence of tumor progression or bile duct obstruction.

¹ Gil-Alzugaray et al. Hepatology, Vol 57, No. 3, 2013.

The typical biochemical picture of REILD is an elevated bilirubin (> 3 mg/dL) in almost all cases, elevated alkaline phosphatase (ALP) and gamma-glutamyl transpeptidase (GGT) in most cases, accompanied by virtually no change in the transaminases (AST and ALT). If liver biopsy is performed, the typical histological appearance is of sinusoidal obstruction that may resemble veno-occlusive disease.

REILD may occur in both non-cirrhotic and cirrhotic patients.

Prophylactic treatment with methylprednisolone and ursodeoxycholic acid starting on the day of SIRT and continued for two months may reduce the incidence of REILD.

In the treatment of REILD, low molecular weight heparin may also be considered but both corticosteroids and heparin may only be useful if commenced very early in the course of the disease. See also *Section 7 Warnings*

10. PATIENT SELECTION AND PRETREATMENT TESTING

Patients with unresectable hepatocellular carcinoma (HCC) or unresectable metastatic liver tumors from primary colorectal cancer refractory to or intolerant of chemotherapy may be considered for treatment with SIR-Spheres. The ability to safely perform resection of tumor(s) is determined by review of individual patient characteristics, for instance during a multi-disciplinary tumor board (MDT). Patients are deemed refractory to or intolerant of chemotherapy after at least one drug or biologic. The potential benefit of selective internal radiation therapy for disease control within the liver is realized in patients with appropriate vascular anatomy of the liver and surrounding tissue, relatively good liver function or functional liver reserve, low lung shunt, and appropriate activity prescription.

10.1 Patient Tests Before Treatment with SIR-Spheres microspheres

The following tests are performed before treatment to ensure the patient is suitable for treatment with SIR-Spheres Y-90 resin microspheres and determine the appropriate dose of SIR-Spheres to prescribe:

- An hepatic angiogram to establish arterial anatomy of the liver
- ^{99m}Tc MAA Scan to determine the percent lung shunting to evaluate lung radiation dose exposure
- Biochemical tests of liver function
- CT, PET/CT, or MRI to determine the extent of the disease, lung, tumor and liver volumes for dosimetry calculations

10.2 Technique for Performing the Intra-Hepatic Technetium MAA Scan

To assess arterial perfusion of the liver and the fraction of radiopharmaceutical tracer that will pass through the liver and lodge in the lungs:

Inject about 150 MBq of ^{99m}Tc MAA into the hepatic artery via a catheter.

Use a large FOV gamma camera and obtain images of the thorax and abdomen (with the same acquisition time).

Draw Region of Interest (ROI) around the whole liver and the whole lung and get the total counts for the lung and the liver.

10.3 Lung Shunt Calculation Procedure

Calculate the lung shunt fraction (L) using following formula:

$$L = \left(\frac{\text{Lung Counts}}{\text{Liver Counts} + \text{Lung Counts}} \right)$$

Equation 1

In order to optimize the risk versus benefit for patients receiving SIR-Spheres microspheres, limiting the radiation exposure to the lungs to ≤ 30 Gy is required. The calculation of estimated radiation exposure to the lungs is given by the following formula:

Activity that may potentially reach the lung:

$$A_{lung} = A_{total} \times L$$

Equation 2

Where:

$$\begin{aligned} A_{lung} &= \text{lung activity [GBq]} \\ A_{total} &= \text{total prescribed activity [GBq]} \\ L &= \text{lung shunt fraction} \end{aligned}$$

The resulting lung dose, given that a given amount of activity shunts from the liver to the lung:

$$D_{lung} = \frac{49670 \times A_{lung}}{M_{lung}}$$

Equation 3

Where:

$$\begin{aligned} D_{lung} &= \text{lung dose [Gy]} \\ A_{lung} &= \text{lung activity [GBq]} \\ M_{lung} &= \text{mass of the lung [g]} \end{aligned}$$

11. CALCULATION OF INDIVIDUAL DOSE

There are two accepted methods for calculating the patient radiation dose, these being the Body Surface Area (BSA) model and the partition model.

11.1 The Body Surface Area (BSA) Method

The BSA method varies yttrium-90 activity according to the size of the patient and the size of the tumour within the liver. The BSA method can be used when treating the volume of a single treated lobe, as well as that of the entire liver.

The approach of lobar treatment versus whole liver treatment with SIR-Spheres microspheres is based on the presence of visible tumors on pre-treatment CT or MR imaging. If liver tumors are only visible in one lobe, then SIR-Spheres microspheres should be administered to that lobe only, thus sparing the contralateral lobe from unnecessary internal radiation.

BSA must first be determined and is calculated from the following equation:

$$BSA = 0.20247 \times H^{0.725} \times W^{0.425}$$

Equation 4

Where:

$$\begin{aligned} BSA &= \text{Body Surface Area} \\ H &= \text{height in meters} \\ W &= \text{weight in kilograms} \end{aligned}$$

11.1.1 BSA prescribed activity calculation for whole liver / bilobar treatment

$$A_{Admin} = (BSA - 0.2) + \left(\frac{V_{Tumor}}{V_{Tumor} + V_{normal\ liver}} \right)$$

Equation 5

Where:

$$\begin{aligned} A_{Admin} &= \text{SIR-Spheres activity to implant [GBq]} \\ V_{Tumor} &= \text{volume of tumor} \\ V_{Normal\ Liver} &= \text{volume of non-tumor tissue in the treated volume} \end{aligned}$$

11.1.2 BSA prescribed activity calculation for lobar or super-selective treatment

In patients who receive lobar or segmental treatment with SIR-Spheres microspheres, the prescribed activity must be reduced in accordance with the size of the portion of the liver being treated.

$$A_{Admin} = \left[BSA - 0.2 + \left(\frac{V_{Tumor}}{V_{Treated}} \right) \right] \times \left[\frac{V_{Treated}}{V_{Liver}} \right]$$

Equation 6

Where:

$$A_{Admin} = \text{SIR-Spheres activity to implant [GBq]}$$

$$\begin{aligned} V_{Tumor} &= \text{volume of tumor in the portion of the liver being treated. (i.e., lobe)} \\ V_{Treated} &= \text{volume of the portion of the liver being treated AND the tumor volume} \\ V_{Liver} &= \text{total volume of the whole liver, including tumor} \\ BSA &= \text{Body Surface Area as per Equation 4} \end{aligned}$$

11.2 Partition Model

This method involves selecting safe radiation doses to the normal liver and lung and implanting the maximum activity that will not exceed these limits. The dose received by the tumor has no upper limit.

The partition model should be used where the tumor mass is a discrete area within the liver. The technique requires two measurements to be made:

1. Measurement of the volume of tumor and normal perfused liver determined from a CT or MR scan
2. Measurement of the proportions of ^{99m}Tc MAA activity that lodges in the tumor, normal perfused liver and lung as determined from imaging.

Therefore, to calculate the activity to be implanted, it is necessary to:

- Determine the volumes of the normal perfused liver and tumor and convert each volume to mass²
- Determine the volume of the lung and convert to mass³
- Using the ^{99m}Tc MAA scan, determine the activity in the lung, tumor and perfused normal liver
- Determine the Tumor to Normal Activity Ratio (TNR) calculated as activity per unit mass of organ or tissue using the following equation:

$$\begin{aligned} TNR &= \frac{A_{Tumor}/M_{Tumor}}{A_{PNL}/M_{PNL}} \\ &= \frac{\text{average counts or activity in tumor}}{\text{average counts or activity in perfused normal liver}} \end{aligned}$$

Equation 7

Where:

$$\begin{aligned} A_{Tumor} &= \text{Activity in tumor} \\ M_{Tumor} &= \text{Mass of tumor} \\ A_{PNL} &= \text{Activity in the perfused normal liver} \\ M_{PNL} &= \text{Mass of the perfused normal liver} \end{aligned}$$

To calculate the total activity to be implanted, use the following equations. The activity required should be calculated to accommodate the lung and normal tissue doses as the limiting factors.

$$A_{Admin} = \frac{D_{PNL} (TNR * M_{Tumor} + M_{PNL})}{49670 * (1 - L)}$$

Equation 8

Where:

$$\begin{aligned} A_{Admin} &= \text{SIR-Spheres activity to implant [GBq]} \\ D_{PNL} &= \text{Mean absorbed dose to perfused normal liver [Gy]} \\ M_{Tumor} &= \text{Mass of tumor [g]} \\ M_{PNL} &= \text{Mass of perfused normal liver [g]} \\ L &= \text{Lung Shunt Fraction} \end{aligned}$$

and,

$$D_{Tumor} = TNR * D_{PNL}$$

Equation 9

Where:

$$D_{Tumor} = \text{Desired absorbed dose to tumor [Gy]}$$

² ICRP, 2015. Radiation Dose to Patients from Radiopharmaceuticals: A Compendium of Current Information Related to Frequently Used Substances. ICRP Publication 128. Ann. ICRP 44(2S).

³ Kappadath SC, Lopez BP, Salem R, Lam MG. Lung shunt and lung dose calculation methods for radioembolization treatment planning. Q J Nucl Med Mol

Imaging. 2021 Mar;65(1):32-42.

12. DOSE PREPARATION PROCEDURE

12.1 Accessories provided by Sirtex:

- V-Vial and capped 21 G needle (SIR-V001)
- V-Vial Holder (SIR-H001)
- Syringe Shield (SIR-S001)

12.2 Additional supplies needed:

- Sterile water for Injection
- Two 25 G vent needles with filters
- 5 mL syringe
- Alcohol swabs
- Forceps
- Dose Calibrator (ion chamber)
- Radiation survey meter
- Radiation spill kit

12.3 Dose Prep Procedure

- Unpack SIR-Spheres microspheres, leaving the glass shipping vial in its lead pot and place on the benchtop.
- Remove the center of the aluminium seal from the sterile v-vial with forceps and wipe the septum with an alcohol swab.
- Place the V-Vial in the V-Vial holder.
- Insert a short 25-gauge needle with filter through the septum of the v-vial until it just pierces the septum to create a vent.
- Leaving the SIR-Spheres microspheres shipping vial in the lead pot and shake to re-suspend the SIR-Spheres microspheres. Resuspension ensures a homogenous solution for dose preparation.
- Open the lead pot and remove the shipping vial using forceps.
- Determine the total activity of SIR-Spheres microspheres in the shipping vial using a dose calibrator, then return the shipping vial to the lead pot.
- Determine the volume to be withdrawn to provide the intended patient specific activity.
- Remove the center of the aluminium seal of the SIR-Spheres microspheres shipping vial, wipe septum with alcohol swab.
- Insert a 25-gauge needle with filter through the septum of the shipping vial to create a vent, ensuring the needle is well clear of the contents in the shipping vial.
- Use a shielded syringe with a 21-gauge needle at least 70 mm long to puncture the septum of the SIR-Spheres microspheres shipping vial, and quickly draw back and forth several times in order to mix the SIR-Spheres microspheres thoroughly.
- Withdraw the precalculated volume.
- Verify the patient dose in the V-Vial by remeasuring the activity in the shipping vial with the dose calibrator, and correct, if necessary.
- Transfer the patient specific activity into the vented V-Vial in V-Vial holder.
- Remove the vent and place the black plug securely into the opening

The patient specific activity is now ready for transport to the SIR-Spheres microspheres implantation room.

13. IMPLANT PROCEDURE

Doctors must refer to the Sirtex Medical Pty Ltd Training Manual for delivering SIR-Spheres resin microspheres before attempting to implant this device.

13.1 Accessories provided by Sirtex:

- Delivery Set (SIR-D001)
- Delivery Box (SIR-B001)
- V-Vial (SIR-V001)
- V-Vial Holder (SIR-H001)

13.2 Additional supplies needed:

- Two 20 mL Luer Lock Syringes filled with non-ionic solution (either 5% glucose/dextrose or water for injection).
 - Do not use saline
 - Patients with diabetes should be monitored for hyperglycemia if 5% dextrose/glucose is used to implant SIR-Spheres
- One 20 mL Luer Lock Syringe filled with non-ionic contrast
- Microcatheter
 - The microcatheters should have an inner diameter of at least 0.021" and a 45° tip configuration for vessels with a 90° angle of origin.
- Sterile absorbent material for protection of angiography suite
- Alcohol swabs
- Forceps or hemostat
- Radiation survey meter
- Radiation spill kit

13.3 Implant Procedure

The hepatic artery catheter is inserted by a trained interventional radiologist under x-ray guidance. This method allows complete control of exactly where the catheter is placed and allows routine checking of the catheter position throughout the implant procedure.

It is essential that SIR-Spheres microspheres are not delivered to other organs, in particular the pancreas, stomach, or duodenum. If there is any possibility of SIR-Spheres microspheres passing down the Gastroduodenal artery (GDA) then the implantation must not proceed. It may be preferable to block shunting vessels with an intraluminal coil or other agent to prevent SIR-Spheres microspheres from flowing to non-target organs.

Note: Virtually all complications from SIR-Spheres microspheres arise from the inadvertent delivery of SIR-Spheres microspheres into small blood vessels that go to the pancreas, stomach or duodenum.

The radiologist must repeatedly check the position of the catheter during the procedure to ensure it remains correctly sited and that reflux of the SIR-Spheres microspheres into other organs does not occur. This is performed by injecting contrast medium through the B line port of the SIR-Spheres Delivery Set during the delivery of SIR-Spheres microspheres. Contrast must not be administered on the D line port.

SIR-Spheres microspheres must be delivered slowly at a rate of no more than 5ml per minute. Rapid delivery may result in a more concentrated suspension of SIR-Spheres that can cause clogging in the microcatheter or reflux back down the hepatic artery and into other organs. At the conclusion of the procedure, the catheter is removed.

14. DISASSEMBLY

- Once the infusion is complete, remove lid from the delivery box.
- Do not disconnect the patient catheter from the tubing set.
- Carefully remove catheter from the patient, while still attached to the Delivery Set. Use care when removing and handling the catheter as it may be contaminated with radioactivity. Coil the catheter and wrap within a sterile towel.
- Leaving everything attached, use a hemostat to push down the needles into the shielded V-Vial and place it (along with all potentially contaminated procedural accessories) into the designated receptacle for radioactive material.
- Radioactive materials must be stored in accordance with local regulations governing radioactive materials storage.

15. CLEANUP AND WASTE DISPOSAL

- Following dose preparation and delivery, survey equipment and people using appropriate radiation survey meter to ensure any contamination is contained.
- If re-usable accessories (i.e., Syringe Shield or Delivery box) are visibly soiled, clean surfaces by spraying with cleaning solution (0.5% sodium hypochlorite, NaOCl) until thoroughly wet. Allow it to stand for 1 minute. Wipe with a soft cloth until all visible soil is removed.
- To disinfect, spray all surfaces with 0.5% NaOCl solution at 6-8 inches until thoroughly wet. Allow to stand for at least 1 minute. Then wipe with a soft cloth to remove any residue.
- Do not use alcohol-based cleaners or abrasives on the reusable accessories.
- Do not attempt to heat sterilize the reusable accessories.
- Materials that are contaminated with radioactive material must be disposed of in accordance with local regulations governing radioactive materials disposal.

16. RADIATION SAFETY

Regulatory and local radiation usage guidelines should be followed concerning implantation and post-implantation care. A properly calibrated radiation survey meter must always be available when handling SIR-Spheres to identify and control potential radioactive contamination risks.

The following are sample measured thermoluminescent dosimetry (TLD) exposures to personnel.

Table 1 – Exposure Dose per Patient for Implant Preparation (Technologist)

	Trunk mSv (mrem)	Lens of the Eye mSv (mrem)	Hands mSv (mrem)
Shallow Dose (0.07mm)	0.027 (2.7)	0.026 (2.6)	0.35 (35)
Deep Dose (10 mm)	0.003 (0.3)	0.004 (0.4)	

Assuming handling of a 3 GBq device and dose preparation time of 30 minutes. TLDs were worn near the pelvis, on the shirt's lapel, and on the working finger.

Table 2 – Exposure Dose per Patient for Implant Procedure (Physician)

	Trunk mSv (mrem)	Lens of the Eye mSv (mrem)	Hands mSv (mrem)
Shallow Dose (0.07mm)	0.038 (3.8)	0.12 (12)	0.32 (32)
Deep Dose (10 mm)	0.004 (0.4)	0.054 (5.4)	

Assuming average patient dose of approximately 2 GBq and dose injection time of 20 minutes.

Exposure data from patients implanted with an average of 2.1 GBq at approximately 5-6 hours post implantation at the following distances from the patient's abdomen:

Table 3 – Post-Implant Exposure

Distance	Exposure
0.25 m	18.8 μ Sv/hr
0.5 m	9.2 μ Sv/hr
1.0 m	1.5 μ Sv/hr
2.0 m	0.4 μ Sv/hr
4.0 m	<0.1 μ Sv/hr

17. CORRECTION FOR DECAY

The physical half-life of yttrium-90 is 64.1 hours. Radioactive decay factors should be applied at the time of patient dose preparation in order to calculate the true value of radioactivity present.

Table 4 – Decay Factors of SIR-Spheres microspheres



Hours	Decay Factor
0.5	0.995
1	0.989
2	0.979
3	0.968
4	0.956
5	0.947
6	0.937
7	0.927
8	0.917
9	0.907
10	0.898
11	0.888
12	0.878
24	0.772

Caution: The time of the initial calibration must be converted to the user's local time.

CH REP

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

SYMBOLS TABLE	
SYMBOL	SYMBOL DEFINITION
	Manufacturer
	Date of manufacture
	Consult instructions for use
	Caution
	Use-by-date
	Lot or batch code
	Catalog number
	Serial number
	Quantity
	Sterilized using irradiation
	Sterilized using steam
	Ionizing radiation
	Single Use Only. Indicates a medical device that is intended for use on a single patient during a single procedure.

SYMBOL	SYMBOL DEFINITION
	Do not resterilize
	Product is not made with natural rubber latex
	Do not use if package is damaged
	Keep dry
	Temperature limit
	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician or licensed healthcare practitioner
	Importer
	Authorized representative in Switzerland
	Authorized representative in the European Community
	CE mark + Notified Body identification number

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