

SIROS™

Control your Y-90 delivery with confidence

Driving Innovation for Y-90 Delivery

SIROS is designed to be simple to use yet versatile enough to allow expanded options for patient-tailored delivery



Visual

- Visually controlled administration enables precision



Intuitive

- Intuitive design for easy set-up, fluid administration and quick disposal



Versatile

- Enhanced control allows for more options to tailor sphere delivery



To see a SIROS demonstration, scan here.

SIROS™

Intuitive Y-90 Delivery



Locking cover

Shields and safely contains the Y-90 during delivery.

Needleless vial

The D-Vial is designed to eliminate touch points and allow more controlled delivery options that may result in a more even distribution of spheres exiting the vial.



Delivery Set

The peel and place tubing set is designed for quick and simple set-up and disposal.



Simple Set-up. As easy as 1-2-3

- 1 Position**
tubing set on base
- 2 Place**
needleless D-Vial holder
- 3 Snap**
on cover to lock and secure SIR-Spheres



For more information,
scan here.

Part Number	Part Name
SIR-10100	SIROS Delivery Dome
SIR-10300	SIROS Delivery Set
SIR-10200	SIROS D-Vial Prep Set

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. **Indications for Use:** SIR-Spheres® Y-90 resin microspheres are indicated for the local tumor control of unresectable hepatocellular carcinoma (HCC) in patients with no macrovascular invasion, Child Pugh-A cirrhosis, well-compensated liver function, and good performance status. They are also indicated for the treatment of unresectable metastatic liver tumors from primary colorectal cancer with adjuvant intra-hepatic artery chemotherapy (IHAC) of FUDR (Floxuridine). **Side Effects:** Common side effects are fever, transient decrease of hemoglobin, mild to moderate abnormality of liver function tests, abdominal pain, nausea, vomiting, and diarrhea. Potential serious effects due to exposure to high radiation include acute pancreatitis, radiation pneumonitis, acute gastritis, acute cholecystitis and radioembolization induced liver disease (REILD). **Consult the Instructions for Use (www.sirtex.com/sir-spheres/risks_adverse-events) for a complete listing of indications, contraindications, side effects, warnings, and precautions.**

Manufacturer

Sirtex Medical Pty Ltd Shop 6, 207 Pacific Highway St Leonards, NSW 2065 Australia Tel: +61 2 9964 8400 Fax: +61 2 9964 8410
Americas Sirtex Medical Inc. 300 Unicorn Park Dr., 2nd FL. Woburn, MA 01801 USA Tel: +1 888 474 7839

www.sirtex.com



@SirtexMedical



company/sirtex-medical-limited



@SirtexMedicalUS

SIR-Spheres® is a registered trademark of Sirtex SIR-Spheres Pty Ltd. ©2026 Sirtex Medical Inc. APM-US-009-05-23 V3

SIRTeX